

	<p><b>Rescue Blanket:</b> Gold and silver used to treat Hyperthermia (heat stroke) OR Hypothermia (cold). Using the silver blanket will keep the casualty warm.</p>		<p><b>Tube gauze bandage:</b> to provide support following soft tissue injury. Available in various sizes. For fingers and limbs.</p>		
	<p><b>Instant disposable Ice Pack:</b> To activate squeeze so that internal liquid bubble breaks. Single use, dispose after use. Useful to reduce pain/swelling/bruising and promote recovery from sprains.</p>		<p><b>Zinc oxide tape:</b> Apply where strong adhesion is required. Do Not apply on sensitive or weak skin.</p>		
	<p><b>Reusable Hot &amp; Cold pack:</b> Cold Use: Reducing swelling, Toothache, Reducing fever. Hot Use; Muscle aches, Back pain Read instructions on package.</p>		<p><b>CPR. Shield / Rescue shield:</b> Only perform CPR if you have been trained to. Dial 111 and wait for instructions.</p>		
	<p><b>Anti-Flamme / Deep Heat</b> Providing deep, penetrating therapy for bumps, bruises, aches and pain. Massage gently on affected area. DO NOT apply on broken skin.</p>	<p><b>For a complete List of Products and Uses visit us on <a href="http://www.nzfirstaid.com">www.nzfirstaid.com</a></b></p> <p><b>4 Main point In First Aid</b></p> <ul style="list-style-type: none"> <li>Start Breathing</li> <li>Stop Bleeding</li> <li>Saturate Burns With Cold Water</li> <li>Support Broken Bones</li> </ul> <p><b>Never attempt to administer more than the basic First Aid if you have not been trained to do so.</b></p>			
	<p><b>Arnica</b> Bruising, sprains and strains, injuries and trauma. Massage gently on affected area. DO NOT apply on broken skin i.e. sores, cuts, blisters.</p>				
	<p><b>Stop Bleed:</b> Ideal for scrapes, grazes and scratches where there is a large surface area and visible bleeding. Also suitable for minor cuts when sprayed onto gauze or tissue and applied to the cut with pressure.</p>				
	<p><b>Optrex:</b> should be used as a general eye wash, useful for removing dust and particles from the eye/s. <b>Refresh:</b> Provide comfortable relief of irritated, dry and tired eyes.</p>				
	<p><b>Hand sanitiser:</b> Apply enough quantity for hand to remain moist for 15 seconds. Rub hand and in between fingers for 15 seconds.</p>	<div style="background-color: #2e8b57; color: white; padding: 10px; border-radius: 10px; text-align: center;"> <p>In an EMERGENCY dial</p> <h1 style="font-size: 48px; margin: 0;">111</h1> <p>Ambulance    Police    Fire</p> <p>Poisons Centre: 0800 764 766</p> </div>			
	<p><b>Gloves:</b> Always wear gloves when attending to a casualty to prevent cross contamination. It is good practice to wear BLUE gloves over plaster when handling food.</p>			<p style="text-align: center;"><b>CUREfast</b> <i>essential First Aid supplies</i></p>	
	<p><b>Splinter Probes:</b> To pry small splinters or glass from the skin. Sterile, single use. Dispose used probes.</p>				